



ITGA2 Antibody, Biotin conjugated

Product Code	CSB-PA011864LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17301
Immunogen	Recombinant Human Integrin alpha-2 protein (188-365AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Integrin alpha-2 (CD49 antigen-like family member B) (Collagen receptor) (Platelet membrane glycoprotein Ia) (GPIa) (VLA-2 subunit alpha) (CD antigen CD49b), ITGA2, CD49B
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ITGA2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.