




ITIH5 Antibody

Product Code	CSB-PA768213LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86UX2
Immunogen	Recombinant Human Inter-alpha-trypsin inhibitor heavy chain H5 protein (35-161AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Inter-alpha-trypsin inhibitor heavy chain H5 (ITI heavy chain H5) (ITI-HC5) (Inter-alpha-inhibitor heavy chain 5), ITIH5, KIAA1953
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	ITIH5

Image



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA768213LA01HU at dilution of 1:100