



ITIH5 Antibody, Biotin conjugated

Product Code	CSB-PA768213LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86UX2
Immunogen	Recombinant Human Inter-alpha-trypsin inhibitor heavy chain H5 protein (35-161AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Inter-alpha-trypsin inhibitor heavy chain H5 (ITI heavy chain H5) (ITI-HC5) (Inter-alpha-inhibitor heavy chain 5), ITIH5, KIAA1953
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	ITIH5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.