



JAK3 Antibody

Product Code	CSB-PA011932LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P52333
Immunogen	Recombinant Human Tyrosine-protein kinase JAK3 protein (351-614AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tyrosine-protein kinase JAK3 (EC 2.7.10.2) (Janus kinase 3) (JAK-3) (Leukocyte janus kinase) (L-JAK), JAK3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	JAK3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.