



KBTBD4 Antibody, Biotin conjugated

Product Code	CSB-PA882119LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NVX7
Immunogen	Recombinant Human Kelch repeat and BTB domain-containing protein 4 protein (151-415AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kelch repeat and BTB domain-containing protein 4 (BTB and kelch domain-containing protein 4), KBTBD4, BKLHD4
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	KBTBD4