



KCNH7 Antibody, HRP conjugated

Product Code	CSB-PA878877LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NS40
Immunogen	Recombinant Human Potassium voltage-gated channel subfamily H member 7 protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Potassium voltage-gated channel subfamily H member 7 (Ether-a-go-go-related gene potassium channel 3) (ERG-3) (Eag-related protein 3) (Ether-a-go-go-related protein 3) (hERG-3) (Voltage-gated potassium channel subunit Kv11.3), KCNH7, ERG3
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	KCNH7