




KCNRG Antibody

Product Code	CSB-PA847638LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N5I3
Immunogen	Recombinant Human Potassium channel regulatory protein (1-229AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Potassium channel regulatory protein (Potassium channel regulator) (Protein CLLD4), KCNRG, CLLD4
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	KCNRG

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA847638LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.