



KCNS2 Antibody, HRP conjugated

Product Code	CSB-PA892470LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9ULS6
Immunogen	Recombinant Human Potassium voltage-gated channel subfamily S member 2 protein (1-183AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Potassium voltage-gated channel subfamily S member 2 (Delayed-rectifier K(+) channel alpha subunit 2) (Voltage-gated potassium channel subunit Kv9.2), KCNS2, KIAA1144
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	KCNS2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.