



KCTD4 Antibody, Biotin conjugated

Product Code	CSB-PA837857LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WVF5
Immunogen	Recombinant Human BTB/POZ domain-containing protein KCTD4 protein (1-259AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	BTB/POZ domain-containing protein KCTD4, KCTD4
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	KCTD4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.