



Keap1 Antibody, HRP conjugated

Product Code	CSB-PA012147LB01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Z2X8
Immunogen	Recombinant Mouse Kelch-like ECH-associated protein 1 protein (1-624AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kelch-like ECH-associated protein 1 (Cytosolic inhibitor of Nrf2) (INrf2), Keap1, Inrf2 Kiaa0132
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Keap1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.