





KIAA0319L Antibody

Product Code	CSB-PA816919LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IZA0
Immunogen	Recombinant Human Dyslexia-associated protein KIAA0319-like protein (479-718AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dyslexia-associated protein KIAA0319-like protein (Adeno-associated virus receptor) (AAVR), KIAA0319L, AAVR KIAA1837
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	KIAA0319L

Image



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA816919LA01HU at dilution of 1:100



Immunofluorescent analysis of A549 cells using CSB-PA816919LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.