



KIAA0319L Antibody, Biotin conjugated

Product Code	CSB-PA816919LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IZA0
Immunogen	Recombinant Human Dyslexia-associated protein KIAA0319-like protein (479-718AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dyslexia-associated protein KIAA0319-like protein (Adeno-associated virus receptor) (AAVR), KIAA0319L, AAVR KIAA1837
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	KIAA0319L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.