



KHDC4 Antibody, Biotin conjugated

Product Code	CSB-PA770364LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z7F0
Immunogen	Recombinant Human KH homology domain-containing protein 4 protein (1-241AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	KH homology domain-containing protein 4 (Brings lots of money 7) (Pre-mRNA splicing factor protein KHDC4), KHDC4, BLOM7 KIAA0907 SNORA80EHG
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	KHDC4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.