



# KHDC4 Antibody, Biotin conjugated

Product Code	CSB-PA770364LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z7F0
Immunogen	Recombinant Human KH homology domain-containing protein 4 protein (1-241AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	KH homology domain-containing protein 4 (Brings lots of money 7) (Pre-mRNA splicing factor protein KHDC4), KHDC4, BLOM7 KIAA0907 SNORA80EHG
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	KHDC4