



KHDC4 Antibody, FITC conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA770364LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q7Z7F0 |
| Immunogen | Recombinant Human KH homology domain-containing protein 4 protein (1-241AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | KH homology domain-containing protein 4 (Brings lots of money 7) (Pre-mRNA splicing factor protein KHDC4), KHDC4, BLOM7 KIAA0907 SNORA80EHG |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Others |
| Target Names | KHDC4 |