



KIAA1257 Antibody, Biotin conjugated

Product Code	CSB-PA891961LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9ULG3
Immunogen	Recombinant Human Uncharacterized protein KIAA1257 protein (1-297AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Uncharacterized protein KIAA1257, KIAA1257
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	KIAA1257
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.