



KIF12 Antibody, Biotin conjugated

Product Code	CSB-PA850294LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96FN5
Immunogen	Recombinant Human Kinesin-like protein KIF12 protein (213-450AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kinesin-like protein KIF12, KIF12
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KIF12
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.