



KLHDC8A Antibody, FITC conjugated

Product Code	CSB-PA816911LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IYD2
Immunogen	Recombinant Human Kelch domain-containing protein 8A protein (1-350AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kelch domain-containing protein 8A (Substitute for delta-EGFR expression 1) (S-delta-E1), KLHDC8A
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	KLHDC8A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.