



KLHL12 Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA700650LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q53G59 |
| Immunogen | Recombinant Human Kelch-like protein 12 protein (1-220AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Kelch-like protein 12 (CUL3-interacting protein 1) (DKIR homolog) (hDKIR), KLHL12, C3IP1 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Others |
| Target Names | KLHL12 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |