




KLHL22 Antibody

Product Code	CSB-PA684469LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q53GT1
Immunogen	Recombinant Human Kelch-like protein 22 protein (2-300AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kelch-like protein 22, KLHL22
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	KLHL22

Image



Western blot
 All lanes: KLHL22 antibody at 8µg/ml
 Lane 1: Jurkat whole cell lysate
 Lane 2: NIH/3T3 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 72, 56 kDa
 Observed band size: 72 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.