



KLHL41 Antibody

Product Code	CSB-PA011995LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60662
Immunogen	Recombinant Human Kelch-like protein 41 protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kelch-like protein 41 (Kel-like protein 23) (Kelch repeat and BTB domain-containing protein 10) (Kelch-related protein 1) (Sarcosin), KLHL41, KBTBD10 KRP1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	KLHL41
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.