



# KLHL41 Antibody

Product Code	CSB-PA011995LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60662
Immunogen	Recombinant Human Kelch-like protein 41 protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kelch-like protein 41 (Kel-like protein 23) (Kelch repeat and BTB domain-containing protein 10) (Kelch-related protein 1) (Sarcosin), KLHL41, KBTBD10 KRP1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	KLHL41