



KLHL42 Antibody, HRP conjugated

Product Code	CSB-PA885803LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P2K6
Immunogen	Recombinant Human Kelch-like protein 42 protein (1-410AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kelch-like protein 42 (Cullin-3-binding protein 9) (Ctb9) (Kelch domain-containing protein 5), KLHL42, KIAA1340 KLHDC5
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	KLHL42
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.