



# KRT33B Antibody, Biotin conjugated

<b>Product Code</b>	CSB-PA621778LD01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q14525
<b>Immunogen</b>	Recombinant Human Keratin, type I cuticular Ha3-II protein (1-404AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Keratin, type I cuticular Ha3-II (Hair keratin, type I Ha3-II) (Keratin-33B) (K33B), KRT33B, HHA3-II HKA3B KRTHA3B
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	KRT33B
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.