



KRT33B Antibody, HRP conjugated

Product Code	CSB-PA621778LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14525
Immunogen	Recombinant Human Keratin, type I cuticular Ha3-II protein (1-404AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Keratin, type I cuticular Ha3-II (Hair keratin, type I Ha3-II) (Keratin-33B) (K33B), KRT33B, HHA3-II HKA3B KRTHA3B
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KRT33B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.