



L2HGDH Antibody, Biotin conjugated

Product Code	CSB-PA864008LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H9P8
Immunogen	Recombinant Human L-2-hydroxyglutarate dehydrogenase, mitochondrial protein (52-463AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	L-2-hydroxyglutarate dehydrogenase, mitochondrial (EC 1.1.99.2) (Duranin), L2HGDH, C14orf160
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	L2HGDH