



LACTB2 Antibody, HRP conjugated

Product Code	CSB-PA706635LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q53H82
Immunogen	Recombinant Human Endoribonuclease LACTB2 protein (6-288AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Endoribonuclease LACTB2 (EC 3.1.27.-) (Beta-lactamase-like protein 2), LACTB2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	LACTB2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.