



LAMA5 Antibody, HRP conjugated

Product Code	CSB-PA012729LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15230
Immunogen	Recombinant Human Laminin subunit alpha-5 protein (3401-3692AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Laminin subunit alpha-5 (Laminin-10 subunit alpha) (Laminin-11 subunit alpha) (Laminin-15 subunit alpha), LAMA5, KIAA0533 KIAA1907
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	LAMA5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.