




LDAH Antibody

Product Code	CSB-PA875685LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H6V9
Immunogen	Recombinant Human Lipid droplet-associated hydrolase protein (1-325AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lipid droplet-associated hydrolase (EC 3.1.1.-) (Lipid droplet-associated serine hydrolase) (hLDAH), LDAH, C2orf43
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	LDAH

Image



Western blot

All lanes: LDAH antibody at 8µg/ml

Lane 1: Hela whole cell lysate

Lane 2: Jurkat whole cell lysate


Lane 3: A431 whole cell lysate

Secondary


Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 38, 33, 23 kDa

Observed band size: 38 kDa



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA875685LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human pancreatic tissue using CSB-PA875685LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.