



LDAH Antibody, Biotin conjugated

Product Code	CSB-PA875685LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H6V9
Immunogen	Recombinant Human Lipid droplet-associated hydrolase protein (1-325AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lipid droplet-associated hydrolase (EC 3.1.1.-) (Lipid droplet-associated serine hydrolase) (hLDAH), LDAH, C2orf43
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	LDAH
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.