



# LDAH Antibody, Biotin conjugated

Product Code	CSB-PA875685LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H6V9
Immunogen	Recombinant Human Lipid droplet-associated hydrolase protein (1-325AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lipid droplet-associated hydrolase (EC 3.1.1.-) (Lipid droplet-associated serine hydrolase) (hLDAH), LDAH, C2orf43
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	LDAH