



CNMD Antibody, Biotin conjugated

Product Code	CSB-PA012854LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75829
Immunogen	Recombinant Human Leukocyte cell-derived chemotaxin 1 protein (215-334AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Leukocyte cell-derived chemotaxin 1 (Chondromodulin) [Cleaved into: Chondrosurfactant protein (CH-SP); Chondromodulin-1 (Chondromodulin-I) (ChM-I)], CNMD, CHMI LECT1 MYETS1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	CNMD
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.