



LHPP Antibody, FITC conjugated

Product Code	CSB-PA867114HC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H008
Immunogen	Recombinant Human Phospholysine phosphohistidine inorganic pyrophosphate phosphatase protein (1-270AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phospholysine phosphohistidine inorganic pyrophosphate phosphatase (hLHPP) (EC 3.1.3.-) (EC 3.6.1.1), LHPP
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	LHPP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.