





LILRA2 Antibody

Product Code	CSB-PA854017LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N149
Immunogen	Recombinant Human Leukocyte immunoglobulin-like receptor subfamily A member 2 protein (9-275AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Leukocyte immunoglobulin-like receptor subfamily A member 2 (CD85 antigen-like family member H) (Immunoglobulin-like transcript 1) (ILT-1) (Leukocyte immunoglobulin-like receptor 7) (LIR-7) (CD antigen CD85h), LILRA2, ILT1 LIR7
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	LILRA2

Image



Western blot
 All lanes: LILRA2 antibody at 4µg/ml + A549 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 53, 52, 47, 49 kDa
 Observed band size: 53 kDa



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA854017LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.