



LILRA4 Antibody, Biotin conjugated

Product Code	CSB-PA012936LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P59901
Immunogen	Recombinant Human Leukocyte immunoglobulin-like receptor subfamily A member 4 protein (289-424AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Leukocyte immunoglobulin-like receptor subfamily A member 4 (CD85 antigen-like family member G) (Immunoglobulin-like transcript 7) (ILT-7) (CD antigen CD85g), LILRA4, ILT7
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	LILRA4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.