



LINC00846 Antibody, Biotin conjugated

Product Code	CSB-PA882112LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NV44
Immunogen	Recombinant Human Putative uncharacterized protein encoded by LINC00846 protein (1-126AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Putative uncharacterized protein encoded by LINC00846, LINC00846, C21orf77
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	LINC00846
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.