



LIX1 Antibody, Biotin conjugated

Product Code	CSB-PA822702LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N485
Immunogen	Recombinant Human Protein limb expression 1 homolog protein (1-282AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein limb expression 1 homolog, LIX1, C5orf11
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.