



LLPH Antibody, FITC conjugated

Product Code	CSB-PA861123LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BRT6
Immunogen	Recombinant Human Protein LLP homolog protein (1-129AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein LLP homolog (Protein LAPS18-like), LLPH, C12orf31 cPERP-G
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	LLPH
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.