



LMO3 Antibody, Biotin conjugated

Product Code	CSB-PA844705LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TAP4
Immunogen	Recombinant Human LIM domain only protein 3 protein (1-145AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	LIM domain only protein 3 (LMO-3) (Neuronal-specific transcription factor DAT1) (Rhombotin-3), LMO3, RBTN3 RBTNL2 RHOM3
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	LMO3