



Loxl1 Antibody

Product Code	CSB-PA013040LA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P97873
Immunogen	Recombinant Mouse Lysyl oxidase homolog 1 protein (95-607AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lysyl oxidase homolog 1 (EC 1.4.3.-) (Lysyl oxidase 2) (Lysyl oxidase-like protein 1), Loxl1, Lox2 Loxl
Immunogen Species	Mus musculus (Mouse)
Research Area	Signal Transduction
Target Names	Loxl1