



Loxl1 Antibody, Biotin conjugated

Product Code	CSB-PA013040LD01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P97873
Immunogen	Recombinant Mouse Lysyl oxidase homolog 1 protein (95-607AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lysyl oxidase homolog 1 (EC 1.4.3.-) (Lysyl oxidase 2) (Lysyl oxidase-like protein 1), Loxl1, Lox2 Loxl
Immunogen Species	Mus musculus (Mouse)
Research Area	Signal Transduction
Target Names	Loxl1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.