



LRRC8E Antibody, Biotin conjugated

Product Code	CSB-PA743982LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6NSJ5
Immunogen	Recombinant Human Volume-regulated anion channel subunit LRRC8E protein (353-630AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Volume-regulated anion channel subunit LRRC8E (Leucine-rich repeat-containing protein 8E), LRRC8E
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	LRRC8E
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.