



LSM6 Antibody

Product Code	CSB-PA013209LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62312
Immunogen	Recombinant Human U6 snRNA-associated Sm-like protein LSm6 protein (1-80AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	U6 snRNA-associated Sm-like protein LSm6, LSM6
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	LSM6