



LSMEM2 Antibody, HRP conjugated

Product Code	CSB-PA822677LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N112
Immunogen	Recombinant Human Leucine-rich single-pass membrane protein 2 protein (1-96AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Leucine-rich single-pass membrane protein 2, LSMEM2, C3orf45
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	LSMEM2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.