



Ly6e Antibody, Biotin conjugated

Product Code	CSB-PA717533LD01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q64253
Immunogen	Recombinant Mouse Lymphocyte antigen 6E protein (21-102AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lymphocyte antigen 6E (Ly-6E) (Stem cell antigen 2) (Thymic shared antigen 1) (TSA-1), Ly6e, Ly67 Sca-2 Tsa-1
Immunogen Species	Mus musculus (Mouse)
Research Area	Cell Biology
Target Names	Ly6e
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.