




LY96 Antibody

Product Code	CSB-PA013254LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y6Y9
Immunogen	Recombinant Human Lymphocyte antigen 96 protein (19-160AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lymphocyte antigen 96 (Ly-96) (ESOP-1) (Protein MD-2), LY96, ESOP1 MD2
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	LY96

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA013254LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.