



MAD2L2 Antibody, HRP conjugated

Product Code	CSB-PA887064LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UI95
Immunogen	Recombinant Human Mitotic spindle assembly checkpoint protein MAD2B protein (1-211AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitotic spindle assembly checkpoint protein MAD2B (Mitotic arrest deficient 2-like protein 2) (MAD2-like protein 2) (REV7 homolog) (hREV7), MAD2L2, MAD2B REV7
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	MAD2L2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.