



MAD2L2 Antibody, HRP conjugated

Product Code	CSB-PA887064LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UI95
Immunogen	Recombinant Human Mitotic spindle assembly checkpoint protein MAD2B protein (1-211AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitotic spindle assembly checkpoint protein MAD2B (Mitotic arrest deficient 2-like protein 2) (MAD2-like protein 2) (REV7 homolog) (hREV7), MAD2L2, MAD2B REV7
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	MAD2L2