



MAGEB18 Antibody, Biotin conjugated

Product Code	CSB-PA853461LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96M61
Immunogen	Recombinant Human Melanoma-associated antigen B18 protein (1-270AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Melanoma-associated antigen B18 (MAGE-B18 antigen), MAGEB18
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	MAGEB18
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.