




MAGEC3 Antibody

Product Code	CSB-PA819466LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TD91
Immunogen	Recombinant Human Melanoma-associated antigen C3 protein (75-197AA)
Raised In	Rabbit
Species Reactivity	Human, Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Melanoma-associated antigen C3 (Cancer/testis antigen 7.2) (CT7.2) (Hepatocellular carcinoma-associated antigen 2) (MAGE-C3 antigen), MAGEC3, HCA2
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer

Image



Western blot

All lanes: MAGEC3 antibody at 2µg/ml

Lane 1: A375 whole cell lysate

Lane 2: Rat thymus tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 72, 39 kDa

Observed band size: 72 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.